

Stirrers for cell cultures and tissue cultures with external Controller 2mag



Product in eShop

Home / [Labortops_Q1_2025](#) / Stirrers for cell cultures and tissue cultures with external Controller 2mag

Product Data

Magnetrührer mit 1/2/3/4 Rührstellen, für zellschonendes Durchmischen von Zellsuspensionen und Kulturbrühen, auch für viskose Kulturen, wartungsfrei, grosse Stellfläche, 100% ruckfrei, gekapseltes Edelstahlgehäuse, wasser-, staub-, keimdicht, IP68, tauchbar, leicht zu reinigen, flache Bauweise, geeignet für Kulturflaschen mit eingehängten, drehenden sowie pendelnden Rührstabsystemen, individuelle Abmessungen auf Anfrage.

Steuergeräte separat bestellen!

Item No.	type no.	stirring posts	speed	Price
17.0114.88	bioMIXdrive1	1	250 U/min	CHF 864.50
17.0114.89	bioMIXdrive2	2	250 U/min	CHF 1'627.60
17.0114.90	bioMIXdrive3	3	250 U/min	CHF 2'366.00
17.0114.91	bioMIXdrive4	4	250 U/min	CHF 3'045.90

Accessories

Item No.	Description	Price
13.1120.00	Magnetic Stirring Bars HUBERLAB.	from CHF 1.40
17.0114.93	Control unit for bioMIXdrive Stirring devices 2mag	from CHF 1'872.00
13.1130.00	Magnetic Stirring Bar Retrievers PTFE Bohlender	from CHF 16.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 24.04.2026 at 05:51 o'clock.