

Microscope BA210E, Motic

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / [Microscope BA210E, Motic](#)

Product Data

Main features

- Infinity Corrected CCIS® Plan Achromatic objectives
- Coaxial coarse and fine focusing
- Quadruple reverse-mounted nosepiece
- Choice between Halogen and LED illumination

- Binocular / Trinocular head Siedentopf type, 30° inclined, 360° rotating (Light split 100:0/20:80)
- Widefield high eyepoint eyepieces N-WF10X/20mm with diopter adjustment on both eyepieces and rubber eyecups
- Reversed quadruple revolving nosepiece
- CCIS EC Plan Achromatic objectives EC PL 4X, 10X, 40X S, 100X S-Oil
- Coaxial coarse and fine focusing system with tension adjustment
- Built in low position coaxial rackless mechanical stage (Right hand control)
- Focusable Abbe condenser N.A. 1.25 with iris diaphragm and slot
- Quartz halogen illumination 6V/30W with intensity control
- 45mm blue filter, immersion oil (5ml), power cord, Allen hexagonal key, thumb screw, vinyl dust cover

Item No.	description	Price
17.7601.01	BA210E Binocular	CHF 1'800.93
17.7601.02	BA210E Trinocular	CHF 2'126.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2026 at 23:02 o'clock.

www.huberlab.ch