

French Square Bottle 60 ml, Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / French Square Bottle 60 ml, Wheaton

Product Data

- Clear USP Type III soda-lime glass
- Ideal for solvent, chemical or sample storage
- Square shape maximizes storage space
- Available with caps attached or bulk packs without caps

Item No.	colour	vol.	d x h	closure size	material	Price
8.2535.51	clear	60 ml	39 x 87 mm	28-400	Cap: White Polypropylene	CHF 110.40
8.2541.19	clear	60 ml	39 x 87 mm	28-400	Cap: White Polypropylene	CHF 153.60
8.2547.73	clear	60 ml	39 x 87 mm	28-400	Cap: Black Phenolic	CHF 160.80
8.2553.32	clear	60 ml	39 x 87 mm	28-400	Cap: Black Phenolic	CHF 132.00
8.2556.68	clear	60 ml	39 x 87 mm	28-400	Cap: White Polypropylene	CHF 575.40
8.2561.12	clear	60 ml	39 x 87 mm	28-400	without cap	CHF 363.84

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 18:31 o'clock.