

Boerner Micro Test Slide Wheaton



Product in eShop

Home / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / Boerner Micro Test Slide Wheaton

Product Data

- Ideal for Boerner-Jones-Lukens microflocculation test
- This slide is used for procedures involving 0.01 to 0.15mL of fluid
- The slide has 10 cells in two rows of five each
- Each cell is 2.5mm deep and is designed to produce flattened drops of uniform depth
- Manufactured from soda-lime glass

Item No.	description	l x w	Price
8.7985.50	Boerner Micro Test Slide	109 x 58 mm	CHF 816.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 07.05.2026 at 03:31 o'clock.