

Lab Mat made of silicone



Product in eShop




Home / [HUBERLAB. green](#) / Lab Mat made of silicone

Product Data

- Made from FDA food grade silicone
- Ideal for every type of laboratory and more specialized areas such as cold rooms, clean rooms, sterile suites and areas where high sterility is a requirement
- Protect countertops from hot items up to 200°C (392°F)
- Chemical resistance: chemically inert and does not react with most chemicals.
- Washable with standard disinfectants in sinks or dishwashers allow for easy clean up for repetitive uses
- Mat can be rolled up for compact storage
- Autoclavable at 121° for 15 minutes

Item No.	colour	dimensions	Price
15.1206.00	front side: solid yellow back side: blue design on grey background	350 x 600 mm	CHF 103.75
15.1206.01	front side: solid purple back side: grey design on purple background	350 x 600 mm	CHF 103.75

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 11.05.2026 at 12:08 o'clock.

www.huberlab.ch