

## Waste Disposal Bag High Transparency, made of PP



 Product in eShop

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Bags](#) / Waste Disposal Bag High Transparency, made of PP

### Product Data

- autoclavable at 121 °C, for autoclaving, the bag must not be completely closed

Item No.	nominal capacity	w x l	film thickness	Price
7.7001.05	2,5 l	200 x 300 mm	0,05 mm	CHF 12.00
7.7001.16	6 l	300 x 500 mm	0,05 mm	CHF 224.80
7.7001.17	22 l	400 x 780 mm	0,05 mm	CHF 258.65
7.7001.25	32 l	600 x 800 mm	0,05 mm	CHF 308.00
7.7001.35	72 l	700 x 1100 mm	0,05 mm	CHF 295.75

### Accessories

Item No.	Description	Price
7.7001.40	Bench-top stand for waste disposal bags ratiolab	from CHF 14.75

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 02.05.2026 at 06:26 oclock.

[www.huberlab.ch](http://www.huberlab.ch)