

Disposable plunging siphon

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Sampling](#) / Disposable plunging siphon

Product Data


- Suitable for use with foodstuffs

Use

- Insert the dip sampler into the liquid and fill it
- Close the thumbhole above the handle using your thumb
- Take out the dip sampler and empty the liquid into the sample bottle

Item No.	length	piercing depth	volume	diameter	Price
15.5387.05	50 cm	35 cm	75 ml	27 mm	CHF 13.10
15.5387.07	75 cm	60 cm	175 ml	27 mm	CHF 24.65
15.5387.10	100 cm	85 cm	275 ml	27 mm	CHF 34.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 29.04.2026 at 09:23 o'clock.