

DESCTR CAB ACRY,12X12X7 EU



Product in eShop

[Home](#) / [Brands](#) / [Sonstige](#) /

Product Data

mit Verschlusslaschen und Scharnier aus Edelstahl, Neopren-Dichtung, zur staub- und feuchtigkeitsfreien Lagerung in Kühl- und Gefrierschränken und auf Labortischen, geeignet für wasserfreie bzw. hygroskopische, biologische und chemische Substanzen sowie elektronische Bauteile und alle Teile, die geregelte Lagerbedingungen erfordern

Item No.	h x w x d	trays	Price
15.5017.01	178 x 305 x 305 mm	2	CHF 679.55
15.5017.02	305 x 305 x 305 mm	3	CHF 885.95
15.5017.03	457 x 305 x 305 mm	4	CHF 1'106.40

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 10:00 o'clock.